## ELECTRONIC APEX LOCATOR EndoEst-Apex 02 (C)

endoline GEO/OFT



www.geosoft.ru/en

## Electronic apex locator EndoEst-Apex 02 (C)



EndoEst-Apex 02 (C) is the most innovative highly accurate electronic apex locator that provides precise measurement of the root canal working length in different chemical conditions due to unique function of data correction which depend on applied intracanal irrigant.

## Features and benefits

High accuracy in the presence of the most popular intra-canal irrigants used in non-surgical endodontic therapy (NaOCl, CHX and EDTA) due to activation of automatic measurement correction.



- Multi-frequency operative system guaranties high accuracy (±0.1 mm) in dry or wet canals.
- A clear Large LCD screen makes it easy to read, and provides a clear display of the file's progressing position along the canal.
- No need for calibration measures accurately in all temperature and moisture conditions inside the canal. It is not necessary to set the device on zero before measuring each individual canal. Self-test function (provided by the apex-locator tester).
- Audible warning system provides a sound-based identification of the file's position along the canal.
- Powered by a rechargeable battery. Auto-power off feature, conserves battery life. Device operating time is approx. 10 hours at one charge.
- Autoclavable file holder to avoid cross infection.
- Special cable fixing protects it from breakage.







## **Accessories**

- File clip "Probe Princh" GE99.121.000;
- Lip clip "Oral Hook" GE99.123.000;
- Cable "Signal Line" (USBB) GE99.169.000.













